



MYCOLIN HOSPITAL DISINFECTANT M 111

MIRIFIC MYCOLIN DISINFECTANT M 111 is a clear soluble phenolic disinfectant but unlike most other is based on a blend of synthetic phenols which act synergistically to achieve an unusually wide spectrum of germicidal activity.

FEATURES

BENEFITS

pH Controlled

MIRIFIC MYCOLIN HOSPITAL DISINFECTANT M 111 has a pH value of 9.5 which, under the normal conditions of dilution falls to around 7.5. The dilutions remain stable over a wide temperature range and retain their germicidal activity over a long period of time.

Investigations over a long period both in the laboratory and of a practical nature have shown MIRIFIC MYCOLIN HOSPITAL DISINFECTANT M 111 to be effective against:

Biocidal Activity

All common Gram positive pathogens.
All common Gram negative pathogens.
Pathogens including *Pseudomonas Aeruginosa*.
Mycobacterium tuberculosis.
Pathogenic fungi.
Viruses of the Pox, Influenza and Herpes.

Safe

MIRIFIC MYCOLIN HOSPITAL DISINFECTANT M 111 has low oral and dermal toxicity.

Directions for use

As a 1.0% Solution – Routine cleaning of walls, floors, operating tables, dressing trolleys and disinfection of sinks, baths sluices etc.

As a 2.0% Solution – Disposal of infected linen, used surgical masks and dressing towels, storage of floor mops, cloth and laboratory brushes. For general disinfection in the post-mortem room.

As a 5.0% Solution – In the pathology laboratory for disposal of highly contaminated materials, cultures etc.

For full toxicological data on MIRIFIC MYCOLIN HOSPITAL DISINFECTANT M 111, please contact your local agent.

Pack Size

Available in 4 x 5 litre cartons.